

Chem 108: Lab Week 13

Sign in / Pick up Papers

Turn in
Global Warming Questions
Scantron

Lab:
Discussion, Experiments & Graded Assignments:

- Week #13 Powerpoint [.html](#), [.ppt](#), Print: [.pdf](#) (6 slides per page)

Experimentation:

- Doing:
Global Warming Questions
Scantron DUE Today

Acid-Base Titration:
Laboratory Manual Acid & Base Titration
(pp. 91-93) Procedures: Part I Completed;
Part II: To Do Today
Part I: Standardization (pg. 94 [DATA form](#))
Part II: Unknown (pg. 95 [DATA form](#), To
Do Today; Complete & have signed
Today)

Complete Report form DUE 15-Nov

Chem 108: Lab

Standardization completed
last week

Individual Titrations of
Unknown
To Do Today



Name: _____
Section: _____

Report Form – Acid Base Titration

Part 1-Standardization of NaOH Solution						
Molarity of HCl used						
Titration	1	2	3	4	5	6
Base buret, final reading (mL)						
Base buret, initial reading (mL)						
Volume of base used (mL)*						
Molarity of NaOH (M)*						
Average molarity of NaOH*					M	

Show the calculations for each of the entries in the Data Table marked with * on the calculations page for one titration.

Neutralization Reactions




There is no difference in the stoichiometry of the reaction of acetic acid & HCl with NaOH, but there is a difference in pH.

Methods for Measuring the pH of an Aqueous Solution



(a) pH paper



(b) Electrodes of a pH meter

Neutralizations / Titrations




Chem 108 titration: phenolphthalein indicator
Chem 120/121 Titration Curves


Aqueous Reactions: Neutralization

pH @ the end point


<http://chemconnections.org/general/movies/pHofSaltSolutions.swf>



The Behavior of Salts in Water			Ion That Reacts with Water
Salt Solution (Examples)	pH	Nature of Ions	
Neutral	7.0	Cation of strong base	None
(NaCl, KBr, Ba(NO ₃) ₂)		Anion of strong acid	
Acidic	<7.0	Cation of weak base	Cation
(NH ₄ Cl, NH ₄ NO ₃ , CH ₃ NH ₃ Br)		Anion of strong acid	
Acidic	<7.0	Small, highly charged cation	Cation
(Al(NO ₃) ₃ , CrCl ₃ , FeBr ₃)		Anion of strong acid	
Basic	>7.0	Cation of strong base	Anion
(CH ₃ COONa, KF, Na ₂ CO ₃)		Anion of weak acid	



<http://chemconnections.org/general/chem108/Acids-Bases%20Guide.html>



<http://chemconnections.org/general/movies/pHEstimation.swf>



pH Estimation

Using the pH benchmarks, estimate the pH of solutions with each of the following hydronium ion concentrations. Click on the boxes to see the estimate and the calculated value of each solution's pH.

pH Benchmarks		
[H ₃ O ⁺]	pH	
1.0 M	0.00	
0.10	1.00	
0.010	2.00	
0.0010	3.00	
1.0 × 10 ⁻⁴	4.00	
1.0 × 10 ⁻⁵	5.00	
1.0 × 10 ⁻⁶	6.00	
1.0 × 10 ⁻⁷	7.00	

[H ₃ O ⁺]	Estimated pH	Calculated pH
0.015 M		
0.058 M		
4.5 × 10 ⁻⁵ M		
3.9 × 10 ⁻² M		

EQUILIBRIUM CO₂ & Oceanic Bicarbonate Buffering



Oceans: pH ~ 8.1 and falling

http://www.tos.org/oceanography/issues/issue_archive/22_4.html

Increasing CO₂ is decreasing ocean pH; long term effects?

http://sos.noaa.gov/datasets/Ocean/ocean_acidification.html


Bicarbonate Buffer Systems



Which is the buffer?

One is bicarbonate buffer to pH = 8 - 9.
One is water plus sodium hydroxide, pH = 8 - 9.

Volunteers to experiment and find out which is which?



Chem 108: Lab

Part 1

Done with a partner.

Name: _____

Section: _____

Report Form - Acid Base Titration

Part 1-Standardization of NaOH Solution



Molarity of HCl used						
Titration	1	2	3	4	5	6
Base buret, final reading (mL)						
Base buret, initial reading (mL)						
Volume of base used (mL)*						
Molarity of NaOH (M)*						
Average molarity of NaOH*					M	

Show the calculations for each of the entries in the Data Table marked with * on the calculations page for one titration.

Completed last week.

Aqueous Reactions: Neutralization

Net Ionic Equations




25.00 mL of $M_{\text{HCL}} = 0.2160 \text{ M}$ nitric acid solution was titrated with a potassium hydroxide solution. It required 24.20 mL as an average of three trials to reach a faint pink color.

$M_{\text{KOH}} = ?$

Aqueous Reactions: Neutralization

Net Ionic Equations



25.00 mL of $M_{\text{HCL}} = 0.2160 \text{ M}$ nitric acid solution was titrated with a potassium hydroxide solution. It required 24.20 mL as an average of three trials to reach a faint pink color.

$$?M_{\text{KOH}} = [M_{\text{HNO}_3} \times V_{\text{HNO}_3} / V_{\text{KOH}}] [? \text{ mol}_{\text{KOH}} / ? \text{ mol}_{\text{HNO}_3}]$$

$$= \frac{0.2160 \text{ mol}_{\text{HNO}_3} \times 0.02500 \text{ L}_{\text{HNO}_3} \times 1 \text{ mol}_{\text{KOH}}}{L_{\text{HNO}_3} \times 0.02420 \text{ L}_{\text{KOH}} \times 1 \text{ mol}_{\text{HNO}_3}} = 0.2231 \text{ M}_{\text{KOH}}$$

QUESTION

A 35.00 mL sample of 0.2250 M HBr was titrated with 42.30 mL of KOH. What is the concentration of the KOH?

- A. 0.0930 M
- B. 0.3030 M
- C. 0.2719 M
- D. 0.1860 M
- E. 0.3720 M


QUESTION

A 35.00 mL sample of 0.2250 M H_2SO_4 was titrated with 42.30 mL of KOH. What is the concentration of the KOH?

- A. 0.0930 M
- B. 0.3030 M
- C. 0.2719 M
- D. 0.1860 M
- E. 0.3720 M

Acid-Base Titration

<https://www.youtube.com/watch?v=9DAB82zLnNE>



	0	2	4	6	8	10	12	14
Methyl violet	Yellow							Violet
Thymol blue	Red		Yellow					Blue
Methyl orange	Red			Yellow				Blue
Methyl red	Red				Yellow			Blue
Bromothymol blue				Yellow		Blue		
Phenolphthalein							Colorless	Pink
Alizarin yellow R							Yellow	Red

Part 2: Unknown Acid

Equipment

From the stockroom:
plastic 1 L bottle
50 mL buret
buret clamp
25 mL vol. pipet and bulb

From the common drawer:
ring stand

From your drawer:
funnel
125 mL flask
250 mL flask
2 beakers (one for waste)
wash bottle

Chem 108: Lab

Part 2: Week 13

To Do (individually) today

Part 2--Determination of Unknown Acid

Unknown code						
Average Molarity of Base from Part 1					M	
Titration	1	2	3	4	5	6
Base buret, final reading (mL)						
Base buret, initial reading (mL)						
Volume of base used (mL)*						
Molarity of unknown acid (M)*						
Average molarity of unknown (M)*					M	

Show the calculations for each of the entries in the Data Table marked with * on the calculations page for one titration.

Take a clean, dry, 125 mL erlenmeyer flask to the stockroom window and get unknown acid solution. Record unknown number. Have data page signed before leaving lab today.

QUESTION

A 35.00 mL sample of hydrochloric acid of unknown concentration was titrated with 42.30 mL of 0.2250 M KOH. What is the concentration of the HCl?

- A. 0.0930 M
- B. 0.3030 M
- C. 0.2719 M
- D. 0.1860 M
- E. 0.3720 M