

## *AA: Quantitative Analysis of Copper Procedure & Data*

### *Preparation and Absorbances of Standard Copper Solutions*

Three standard solutions are needed for a standard curve. A standard stock solution, which contains 5.0 mg of  $\text{Cu}^{2+}$  per liter of solution, has been prepared. You are to prepare two additional concentrations by serial dilution.: 2.50 ppm, and 1.25 ppm. An additional point on the standard curve is the origin, 0.0 mg/liter, which is the blank.

Check out two 50-mL volumetric flasks, two 25.00 mL pipets and a pipet bulb. Prepare the two solutions of 2.50 ppm and 1.25 ppm respectively. Save them.

Each group will develop a calibration curve from the absorbance of each of the standard concentrations.

---

### *Individual Sample Preparation and Analysis*

The copper containing tablet will be treated so that all of the copper is in solution and a serial dilution made to produce a solution with a low enough  $\text{Cu}^{2+}$  concentration to be read in the atomic absorption spectrophotometer.

Check out a mortar & pestle, two 100-mL volumetric flasks, two 10.00 mL pipets and a pipet bulb.

1. Dissolve the copper in the tablet.

Grind your unknown tablet in the mortar and pestle. Transfer to a 100 mL beaker. Add ~ 70 mL of 3 N  $\text{HNO}_3$  and stir for ~5 minutes. Gravity filter directly into a 100.00 mL volumetric flask. and dilute to the mark with deionized water.

2. Dilute the tablet's solution.

Using a 10.00 mL pipet, place a 10.00 mL aliquot of the tablet solution into a second 100.00 mL volumetric flask. Add deionized water to dilute to 100 mL.

3. Determine the absorbance of the diluted vitamin solution. Record on the Report Form.